NIEUWE VARIANTEN **NORA LUNATEC FUSION**

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nora® LUNATEC FUSION

► Application example 1

Soft cushioning insole (perforated)

Lunatec fusion 30 perforated

Thickness: 4 mm

Function: Base and cushioning layer for long sole insole

Lunatec fusion 40 smooth

Thickness: 8 mm

Function: Stabilisation in rearfoot area

1. Roughen and heat

Roughen both materials in contact areas, sand rearfoot stabilisation to wedge shape and then **heat** in oven at 130°C//266°F.

Heating times:

8 mm approx. 4 minutes
4 mm approx. 1.5 minutes


2. Deep draw and allow to cool

Place materials on last, deep draw and make sure that sufficient pressure is applied to the drawing bladder. Press material firmly during the deep drawing process to prevent air insertions.

Cooling time: approx. 8 minutes

3. Grinding

Grind insoles to shape.

i

Tip: Depending on the treatment, a thermoplastic reinforcing material (EVA-based) can also be inserted between the layers, e.g. for heavier body weights or higher loading.

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
► Application example 2

Soft cushioning insole

	Lunatec fusion 30 smooth Thickness: 6 mm Function: Forefoot cushioning
	Lunatec fusion 40 smooth Thickness: 8 mm Function: Stabilisation in rearfoot area
	Lunatec fusion 40 perforated Thickness: 8 mm Function: Stabilisation with high arch

1. Grind and heat

Grind **Lunatec fusion 30** well in contact area and grind **Lunatec fusion 40** (smooth) in wedge shape to avoid edges and then heat in oven at 130°C//266°F.

Heating times:

8 mm	approx. 4 minutes
(slightly less with perforated material)	
6 mm	approx. 3 minutes


2. Deep draw and allow to cool

Place all materials correctly on last, deep draw and make sure that sufficient pressure is applied to the drawing bladder. Press material firmly during the deep drawing process to prevent air insertions.

Cooling time: approx. 8 minutes

3. Grinding

Grind insoles to shape.

i **Tip:** With this type of insole manufacture, it is important to cleanly sand the contact surfaces of the half sole wedge-shaped to avoid edges and ensure ideal bonding.

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Lunatec fusion 30 smooth

Thickness: 4 mm

Function: Base and cushioning layer for long sole bedding

Lunatec fusion 40 perforated

Thickness: 8 mm

Function: Stabilisation in rearfoot area

1. Roughen and heat

Roughen both materials in contact areas, grind rearfoot stabilisation to wedge shape and then **heat** in oven at 130°C // 266°F.

Heating times:

8 mm just under 4 minutes
4 mm approx. 2 minutes

2. Deep draw and allow to cool

Place materials on last, deep draw and make sure that sufficient pressure is applied to the drawing bladder. Press material firmly during the deep drawing process to prevent air insertions.

Cooling time: approx. 8 minutes

3. Grinding

Grind insoles to shape.

Tip: Depending on the treatment, a thermoplastic reinforcing material (EVA-based) can also be inserted between the layers, e.g. for heavier body weights or higher loading.

6

Lunatec fusion 30
378 blauw

* C 1 3 3 1 4 1 2 *

Lunatec fusion 40
09 wit


* C 1 3 3 1 4 7 8 *

Lunatec fusion 20
60 lichtgrijs

* C 1 3 3 1 3 5 2 *

Lunatec fusion 50
27 lichtblauw

* C 1 3 3 1 5 0 8 *